

CSM/CVSM SERIES DEWARS

SUPERCONDUCTING MAGNET DEWARS


CSM

U171006


CSM LN₂ SHIELDED DEWARS




Model	Inside Diameter (Inches)	Outside Diameter (Inches)	Inner Depth (Inches)	Overall Height (Inches) ₁
CSM648	6.00	12.00	48.00	54.00
CSM7548	7.50	13.00	48.00	54.00
CSM8548	8.50	14.00	48.00	54.00
CSM9548	9.50	16.00	48.00	54.00
CSM10548	10.50	16.00	48.00	54.00
CSM12548	12.50	18.00	48.00	54.00
CSM1548	15.00	20.00	48.00	54.00



MAGNET DEWARS



MAGNET DEWAR, LN₂ SHIELDED WITH BELLY



MAGNET DEWAR, LN₂ SHIELDED WITH BELLY AND TAIL